

## How to write a PaR-PaR script

PaR-PaR is a high-level language that enables Biologists to more quickly design experiments that utilize liquid-handling robotics.

General guidelines for PaR-PaR scripts:

- All lines beginning with the pound character ‘#’ are ignored by PaR-PaR. These lines can contain comments that help explain/document what is being performed in the script.
- Lines may be either tab- or space-separated. For this reason, the separated elements within a line (such as a plate name alias) should not themselves contain tab or space characters (*i.e.*, use “DrinksPlate” rather than “Drinks Plate”).
- Variable definitions (*e.g.*, plate name aliases) must precede their use in the script.

A PaR-PaR configuration file consists of:

1. A link to robotic table file. In most instances, table files are created using software distributed with the robot. There may also be a set of ready-made table files to choose from.
2. The script itself, consisting of two logical sections: definitions and actions.

### Definitions Section

#### (Experiment) name

The name of the experiment may optionally be specified:


```
NAME      BreakfastDrinks
```

#### Table (file)

The name of the robotic table file to use must be specified.

```
TABLE      BreakfastDrinks.ewt
```

For the stand-alone version of PaR-PaR, it is necessary to include the (.ewt) table file in the “Tables” folder inside the PaR-PaR folder:



It is very important to verify that the correct table file is specified for the experiment.

## Documentation Section

A documentation section, enclosed by a pair of tripled-quotation mark characters ("\""), may be included in the script. This section is operationally ignored by PaR-PaR (much like lines beginning with the pound character '#'), but can help explain/document what is being performed in the script.

Here is an example documentation section:

```
"""
Recipe for breakfast drinks.
"""
```

## Plate (name aliases)

Aliases may be specified for the plate names (that are themselves specified in the robotic table file):

#	alias	name
PLATE	DrinksPlate	PL4

Note: plate name aliases must not contain any space characters.


## Locations (sources and destinations)

Locations specify plates and wells. Wells are specified as follows:

1. For a single well: by letter-number notation (*e.g.*, "A1") or by number-only notation (*e.g.*, "1").
2. For multiple consecutive wells: the first well, then the character '+', and then the total number of wells. For example, "A1+4" indicates the four consecutive wells starting from well A1 (*i.e.*, wells A1, B1, C1 and D1).
3. For multiple non-consecutive wells: the wells separated by commas ',' (*e.g.*, "A1, C1, E1, G1").
4. Consecutive and non-consecutive wells may be interspersed (*e.g.*, "A1+4, F1").


A plate and its wells are separated by a colon ':' (*e.g.*, "PL8:A1+4, F1"). Multiple plates are separated by forward slashes '/' (*e.g.*, "PL8:A1+4, F1/PL5:A1, C1, E1"). Plate name aliases may be used instead of the plate name (*e.g.* "WaterPlate:A1+4, F1" instead of "PL8:A1+4, F1").

Note: locations must not contain any space characters.




**Figure 1.** Plate well locations in letter-number notation (left), and in number-only notation (right).

Wells are ordered sequentially first from top to bottom, and then left to right. 'A1+4' (or '1+4'), indicates wells 'A1, B1, C1, D1' (or '1, 2, 3, 4').



**Figure 2.** Plate well locations indicated by 'A1+4' (left) or '1+4' (right).



**Figure 3.** Plate well locations indicated by 'B1+4' (letter-number notation, left) or '2+4' (number-only notation, right). Note that for this particular plate geometry (4 rows by 6 columns), well location A2 (or 5) follows well location D1 (or 4).

## Methods

Methods specify the type (or class) of a liquid, as well as the pipetting method for transferring the liquid from one location to another.

The method should be one of the following:

1. LC\_W\_Bot\_Bot      Water, aspirated from the bottom of the well, dispensed to the bottom.
2. LC\_W\_Lev\_Bot      Water, aspirated from the liquid meniscus level, dispensed to the bottom.
3. LC\_W\_Lev\_Lev      Water, aspirated from the liquid level, dispensed to the liquid level.
4. LC\_W\_Lev\_Air      Water, aspirated from the liquid level, dispensed in air above the liquid.
5. DEFAULT            LC\_W\_Bot\_Bot; see immediately below.

Outside of COMPONENT definition statements (see immediately below), DEFAULT refers to the method specified for the component in its definition statement. If there is no method pre-specified for the component, the method defaults to LC\_W\_Bot\_Bot.

## Component(s)

Components may be specified with names, locations, and pipetting methods:

#	name	location	method
COMPONENT	Water	PL8:A1+4, F1	LC_W_Lev_Air

PL8      Water

	1	2	3	4	5	6	7	8	9	10	11	12
A	A1	A2	A3	A4	A5	A6	A7	A8	A9	A10	A11	A12
B	B1	B2	B3	B4	B5	B6	B7	B8	B9	B10	B11	B12
C	C1	C2	C3	C4	C5	C6	C7	C8	C9	C10	C11	C12
D	D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	D11	D12
E	E1	E2	E3	E4	E5	E6	E7	E8	E9	E10	E11	E12
F	F1	F2	F3	F4	F5	F6	F7	F8	F9	F10	F11	F12
G	G1	G2	G3	G4	G5	G6	G7	G8	G9	G10	G11	G12
H	H1	H2	H3	H4	H5	H6	H7	H8	H9	H10	H11	H12

**Figure 4.** Component location PL8 : A1+4, F1 on the sample plate (8 rows by 12 columns).

Note: component names must not contain any space characters.

## Volume (aliases)

Aliases may be specified for volumes (in  $\mu\text{L}$ ):

#	alias	volume ( $\mu\text{L}$ )
VOLUME	DrinkVol	50
VOLUME	WaterVol	25

Note: volume aliases must not contain any space characters.

## Recipe(s)

Recipes may be defined to alias specific mixtures of components, such as a PCR master mix.

The definition of a recipe begins with a `RECIPE` line that specifies the name of the recipe as a whole. The lines that immediately follow specify the names and contents of the sub-recipe(s) associated with the recipe (one sub-recipe per line). In the example below, the recipe `Drinks` has three associated sub-recipes: `chai`, `coffee` and `lemonade`. Each sub-recipe name is followed by a colon ':', which in turn is followed by the contents (specified in pairs of component name or locations, and volume (in  $\mu\text{L}$ )) of the sub-recipe. Plate name aliases, components, and volume aliases (after they have been defined) may be used when defining a recipe.

#	name					
RECIPE	Drinks					
#	component1	volume1	component2	volume2	component3	volume3
chai:	TeaExtract	30	Syrup	30	Water	WaterVol
coffee:	BeanExtract	30	Milk	30	Water	WaterVol
lemonade:	LemonJuice	15	PL7:18	45	Water	WaterVol

## Actions Section

### Options

Options specify what should happen after liquid has been dispensed into a well.

Each option should be one of the following:

1. **MIX:volume×repetitions** Aspirate/dispense *volume* µL *repetitions* times (e.g. MIX:25×20)

Note: options must not contain any space characters.

### Make (a recipe)

The MAKE action prepares a defined recipe, or sub-recipe(s), at the specified location(s).

Sub-recipes are prepared in separate locations, and thus a location must be specified for each sub-recipe. For consecutive well locations, the sub-recipes are prepared sequentially into consecutive wells.

#	recipe:sub-recipe	location	method	options
MAKE	Drinks	DrinksPlate:A6+3	DEFAULT	MIX:25×20
MAKE	Drinks:coffee,lemonade	DrinksPlate:A1+2	DEFAULT	MIX:30×10

Note: recipe:sub-recipe(s) must not contain any space characters.

### Spread (a component)

The SPREAD action distributes a single defined component (or the same liquid present in one or more source locations) to one or more destinations.

#	component	destination	volume (uL)	method	options
SPREAD	Water	PL6:A4+10,A6	DrinkVol	DEFAULT	MIX:25×20
SPREAD	PL4:A1+4	PL6:A4+10,A6	DrinkVol	DEFAULT	MIX:25×20

### Transfer (liquids)

The TRANSFER action distributes liquids one-to-one from source to destination locations. A destination location must be specified for each source location.

#	source	destination	volume (uL)	method	options
TRANSFER	PL1:A1+3	PL6:A7+3	150	LC_W_Bot_Bot	MIX:15×8

## Putting it all together: an example PaR-PaR script

```
NAME      BreakfastDrinks
TABLE     BreakfastDrinks.ewt
```

```
""
Recipe for breakfast drinks.
""
```

```
#      alias      name
PLATE  DrinksPlate PL4
```

```
#      name      location      method
COMPONENT Water      PL8:A1+4,F1  LC_W_Lev_Air
COMPONENT TeaExtract  PL7:17      LC_W_Lev_Bot
COMPONENT Syrup       PL7:18      LC_W_Lev_Bot
COMPONENT Milk        PL7:19      LC_W_Lev_Bot
COMPONENT BeanExtract  PL7:20      LC_W_Lev_Bot
COMPONENT LemonJuice   PL7:21      LC_W_Lev_Bot
```

```
#      alias      volume(uL)
VOLUME DrinkVol    50
VOLUME WaterVol     25
```

```
#      name
RECIPE Drinks
#      component1  volume1  component2  volume2  component3  volume3
chai:   TeaExtract  30        Syrup       30        Water       WaterVol
coffee: BeanExtract 30        Milk        30        Water       WaterVol
lemonade: LemonJuice 15        PL7:18      45        Water       WaterVol
```

```
#      recipe:sub-recipe      location      method      options
MAKE  Drinks                  DrinksPlate:A6+3  DEFAULT     MIX:25x20
MAKE  Drinks:coffee,lemonade  DrinksPlate:A1+2  DEFAULT     MIX:30x10
```

```
#      component  destination  volume(uL)  method      options
SPREAD  Water      PL6:A4+10,A6  DrinkVol    DEFAULT     MIX:25x20
```

```
#      source      destination  volume(uL)  method      options
TRANSFER PL1:A1+3      PL6:A7+3    150         LC_W_Bot_Bot  MIX:15x8
```